

# Magnetic 3D Bioprinting: A novel high-throughput and high-content assay for toxicity screening

Hubert Tseng<sup>1</sup>, Jacob A. Gage<sup>1</sup>, William L. Haisler<sup>2</sup>, Shane K. Neeley<sup>2</sup>, Glauco R. Souza<sup>1</sup>, T. C. Killian<sup>3</sup>

<sup>1</sup>Nano3D Biosciences, Houston, TX, <sup>2</sup>Department of Bioengineering, Rice University, Houston, TX, <sup>3</sup>Department of Physics, Rice University, Houston, TX



## Background

A growing demand exists for three-dimensional (3D) *in vitro* assays that mimic the tissue of interest to predict the human *in vivo* response with high throughput. Animal models are representative of the native tissue but expensive and low-throughput. Currently available *in vitro* assays are rapid but poor representatives of the native tissues, as they are typically two-dimensional (2D) models on rigid substrates. **Thus, the choice of assay becomes a tradeoff between efficiency and accuracy, leaving an unmet need for an assay system that is both representative and high-throughput.**


To that end, we introduce magnetic 3D bioprinting for high-throughput screening. The basis is the magnetization of cells by binding magnetic nanoparticles to them. After resuspension in media, these cells can be rapidly printed into particular shapes, like rings or dots/spheroids, with the use of magnetic forces in high-throughput formats. These printed structures immediately demonstrate a dose-dependent response, which can be visually monitored. An iPod-based system is used for imaging, which is programmed to image whole plates at specific intervals, thereby forgoing the need to efficiently image well-by-well under a microscope.

In this study, we have applied magnetic 3D bioprinting to generate a basic toxicity model using 3T3 murine embryonic fibroblasts, and a specific model for vascular toxicity using A10 vascular smooth muscle cells.

**Hypothesis: Magnetic 3D bioprinting can be used to rapidly print cells into structures that mimic native tissue for high-throughput compound screening**

## Magnetic 3D Bioprinting

- Cells are incubated with NanoShuttle (Nano3D Biosciences) overnight
- The next day, cells are levitated to induce synthesize ECM for a few hours
- Cultures are then broken apart and printed into rings for 15 min - 1 h (1.50K cells/ring, 75K cells/dot) in 96-well plates →
- Magnetic field removed and cells are allowed to close



**Magnetic 3D bioprinting can rapidly and simultaneously print multiple tissue-like structures**


## iPod-Based Imaging System




- Images of 3D printed tissues are taken with an iPod (Apple Computer)
- iPod is programmed using a freely available app (Experimental Assistant, Nano3D Biosciences) to image in real-time
- ← Imaging setup fits within a standard incubator
- iPod imaging forgoes time-consuming well-by-well imaging with a microscope

**Imaging with an iPod improves throughput and efficiency**




## Results



Contraction and dilation in rings and dots after printing and removing the magnetic field




Contraction of dots of 3T3 murine embryonic fibroblasts exposed to retinoic acid over 5 h. As retinoic acid concentration increases, the dots are unable to shrink as much. Scale bar = 5 mm.



**The response of magnetically 3D bioprinted models correlates with viability and cytoskeletal organization**

## Results

### In Vitro Aortic "Ring" Assay



Results of the aortic ring assay with A10 rat vascular smooth muscle cells exposed to blebbistatin. **(A)** A10 rings before and 5 and 300 min after adding blebbistatin to the rings. Scale bar = 5 mm. **(B)** Ring contraction/dilation over 300 min. **(C)** Area change after 5 min as a function of concentration. Solid line = 0, dotted line = control. \* $p < 0.05$  v. control. **(D)** IHC stains for phosphorylated myosin light chain (pMLC, red) 5 min after adding blebbistatin. Black scale bar = 5 mm, white scale bar = 50  $\mu\text{m}$ . Note with increasing concentrations of blebbistatin, the rings dilate and don't stain for pMLC, which is associated with smooth muscle contraction.<sup>3</sup>

**Magnetic 3D bioprinting can be used to create organ-specific models for high-throughput screening**

## Conclusions

- Magnetic 3D bioprinting rapidly prints multiple tissue-like structures that respond to compounds within 24 h
- Ring/dot shrinkage correlates with viability, cytoskeletal organization, and cell behavior
- Magnetic 3D bioprinting can be applied specifically to certain systems, like the cardiovascular system, to create organotypic models for high-throughput screening
- iPod-based system improves throughput and efficiency

## References

1. Haisler, W. L. et al. Three-dimensional cell culturing by magnetic levitation. *Nat. Protoc.* 8, 1940-9 (2013).
  2. Timm, D. M. et al. A high-throughput three-dimensional cell migration assay for toxicity screening with mobile device-based macroscopic image analysis. *Sci. Rep.* 3, 3000 (2013). →
  3. Tan, J. L. et al. Control of nonmuscle myosins by phosphorylation. *Annu. Rev. Biochem.* 61, 721-59 (1992).
- Scan the QR-codes in the poster for videos of the printed structures and our publication in *Scientific Reports!*

For more information, go to [www.n3dbio.com](http://www.n3dbio.com), or email us at [info@n3dbio.com](mailto:info@n3dbio.com)

## Acknowledgements/Funding

This work was supported by the National Science Foundation [Small Business Innovation Research Phase I 0945954 and Phase II 1127551]; and the State of Texas Emerging Technology Fund. The authors would like to thank Ahmed Nasser Salameh, Ph.D., Alexes Daquinag, Ph.D., and Mikhail Kolonin, Ph.D., University of Texas Health Science Center at Houston, for donating the 3T3s