

# A novel vascular "ring" assay for smooth muscle contractility using magnetic 3D bioprinting



Hubert Tseng<sup>1</sup>, Jacob A. Gage<sup>1</sup>, William L. Haisler<sup>1</sup>, Shane K. Neeley<sup>1</sup>, T. C. Killian<sup>2</sup>, Herbert G. Barthlow<sup>3</sup>, Matthew Wagoner<sup>3</sup>, Glauco R. Souza<sup>1</sup>  
<sup>1</sup>Nano3D Biosciences, Houston, TX, <sup>2</sup>Department of Physics, Rice University, Houston, TX, <sup>3</sup>AstraZeneca, Waltham, MA

## Background

One of the most common assays for studying pharmacodynamics on vascular smooth muscle is wire myography, in which *ex vivo* vessel segments are exposed to a drug and its contraction and contractile forces are recorded. This assay suffers from high costs and low throughput, limiting its widespread use. Yet *in vitro* alternatives are insufficient, as they are commonly performed on stiff two-dimensional (2D) substrates that poorly mimic the soft vasculature. **Thus, there is an unmet need for an *in vitro* assay that mimics the 3D structure of the blood vessel and is high-throughput.**


To address this unmet need, we apply magnetic 3D bioprinting towards assaying vascular toxicity for high-throughput screening as an *in vitro* alternative to wire myography. The basis is the magnetization of cells by binding magnetic nanoparticles to them.<sup>1</sup> These cells can be rapidly printed into rings using magnetic forces in high-throughput formats.<sup>2</sup> These rings visually mimic vessel segments and immediately respond to compounds by contracting or dilating. Contraction is captured using an iPod-based system programmed to image whole plates at specific intervals, thereby forgoing the need to image well-by-well under a microscope.

In this study, we have validated this vascular "ring" assay using A10 rat vascular smooth muscle cells and primary aortic smooth muscle cells (ASMC).

**Hypothesis: Magnetic 3D bioprinting can be used to rapidly print cells into vascular "rings" that contract or dilate *in vitro* for high-throughput testing.**

## Magnetic 3D Bioprinting

- Cells are incubated with NanoShuttle (Nano3D Biosciences) overnight
- The next day, cells are levitated to induce synthesis ECM for a few hours
- Cultures are then broken apart and printed into rings for 6 h (200K cells/ring) in 96-well plates →
- Magnetic field removed and cells are allowed to contract



**Magnetic 3D bioprinting can rapidly and simultaneously print multiple vascular rings**


## iPod-Based Imaging System




- Images of 3D printed rings are taken with an iPod (Apple Computer)
- iPod is programmed using a freely available app (Experimental Assistant, Nano3D Biosciences) to image in real-time
- ← Imaging setup fits within a standard incubator
- iPod imaging forgoes time-consuming well-by-well imaging with a microscope

**Imaging with an iPod improves throughput and efficiency**




## Results



Schematic of contraction and dilation of vascular rings immediately after printing





Contraction of A10 rat vascular smooth muscle cells exposed to blebbistatin. (A) A10 rings before and 5 and 300 min after adding blebbistatin to the rings. Scale bar = 5 mm. (B) Ring contraction/dilation as a function of time.



Graphs: Solid line indicates zero. Dotted line indicates control value. \*:  $p < 0.05$  v. control. IHC: Blue = DAPI, Red = pMLC, Black scale bar = 5 mm. White scale bar = 50 µm.

## Results



Dose-response of A10s and ASMCs exposed to norepinephrine, as measured by area change after 5 min (A, C) and 5 h (D), and (B, E) pMLC staining. A10s contracted more with higher concentrations, while ASMCs contracted in a non-monotonic manner over a longer time.

## Conclusions

- Magnetic 3D bioprinting rapidly prints multiple vessel-like rings that contract/dilate to compounds within 24 h
  - Contraction is a label-free metric where rings are available for post-assay experimentation or high-content imaging
  - iPod-based system improves throughput and efficiency
  - Results show differences in response between A10s and ASMCs
- Magnetic 3D bioprinting can be applied towards high-throughput vascular toxicity screening**

## References

1. Haisler, W. L. et al. Three-dimensional cell culturing by magnetic levitation. *Nat. Protoc.* 8, 1940-9 (2013).
2. Timm, D. M. et al. A high-throughput three-dimensional cell migration assay for toxicity screening with mobile device-based macroscopic image analysis. *Sci. Rep.* 3, 3000 (2013). →

Scan the QR-codes in the poster for videos of the printed structures and our publication in Scientific Reports!

For more information, go to [www.n3dbio.com](http://www.n3dbio.com), or email us at [info@n3dbio.com](mailto:info@n3dbio.com)

## Acknowledgements/Funding

This work was supported by the National Science Foundation [Small Business Innovation Research Phase I 0945954 and Phase II 1127551]; and the State of Texas Emerging Technology Fund.

