

Magnetic 3D Bioprinting: A novel high-throughput and high-content assay for toxicity screening



Hubert Tseng¹, Jacob A. Gage¹, William L. Haisler², Shane K. Neeley², Glauco R. Souza¹, T. C. Killian³
¹Nano3D Biosciences, Houston, TX, ²Department of Bioengineering, Rice University, Houston, TX, ³Department of Physics, Rice University, Houston, TX

Background

A growing demand exists for three-dimensional (3D) *in vitro* assays for toxicity screening. Animal models are representative of the native tissue but expensive, low-throughput, and ultimately, not human. Currently available *in vitro* assays are rapid but poorly representative tissues, as they are typically two-dimensional (2D) models on rigid substrates. **Thus, the choice of assay becomes a tradeoff between efficiency and accuracy, leaving an unmet need for an assay system that is both representative and high-throughput.**


To that end, we introduce magnetic 3D bioprinting for high-throughput screening. The basis of this approach is the magnetization of cells by binding magnetic nanoparticles to them.¹ After resuspension in media, these cells can be rapidly and simultaneously printed into particular shapes, like rings or spheroids, with the use of magnetic forces.² These printed structures immediately demonstrate a dose-dependent response, which can be visually monitored. An iPod-based system is used for imaging, which is programmed to image whole plates at specific intervals, thereby forgoing the need to efficiently image well-by-well under a microscope.

In this study, we have applied magnetic 3D bioprinting to generate a basic toxicity model using 3T3 murine embryonic fibroblasts, and a specific model for vascular toxicity using vascular smooth muscle cells.

Hypothesis: Magnetic 3D bioprinting can be used to rapidly print cells into structures that mimic native tissue for high-throughput compound screening


Magnetic 3D Bioprinting

- Cells are incubated with NanoShuttle (Nano3D Biosciences) overnight
- The next day, cells are levitated to induce ECM synthesis for a few hours
- Cultures are then broken apart and printed into rings for 15 min - 6 h (150K cells/ring, 75K cells/spheroid) in 96-well plates →
- Magnetic field removed and cells are allowed to close



Magnetic 3D bioprinting can rapidly and simultaneously print multiple tissue-like structures


iPod-Based Imaging System




- Images of 3D printed tissues are taken with an iPod (Apple Computer)
- iPod is programmed using a freely available app (Experimental Assistant, Nano3D Biosciences) to image in real-time
- ← Imaging setup fits within a standard incubator
- iPod imaging forgoes time-consuming well-by-well imaging with a microscope

Imaging with an iPod improves throughput and efficiency


Results




Contraction and dilation in rings and dots after printing and removing the magnetic field



Contraction of spheroids of 3T3 murine embryonic fibroblasts exposed to retinoic acid over 5 h. As retinoic acid concentration increases, the dots are unable to contract as much. Scale bar = 5 mm.



(A) Contraction of 3T3 spheroids exposed to over 5 h and (B) the rate of contraction over 3 h as a function of concentration. Spheroid contraction demonstrates dose-dependence.






Viability of contracted spheroids exposed to retinoic acid, doxorubicin, and forskolin. As expected, at higher concentrations, the spheroids are less viable, matching spheroid size and the ability to contract. Scale bar = 500 μm

Magnetically 3D bioprinted models measure toxicity in 3D environments, which correlates with viability

Results

In Vitro Aortic "Ring" Assay



Conclusions

- Magnetic 3D bioprinting rapidly prints 3D structures that immediately respond to compounds in a manner related to viability
- High-content experimentation can be performed to explore mechanisms of action and yield more data.
- Organotypic models, like of vasculature, can be created for high-throughput screening
- iPod-based system improves throughput and efficiency

References

1. Haisler, W. L. et al. Three-dimensional cell culturing by magnetic levitation. Nat. Protoc. 8, 1940-9 (2013).
2. Timm, D. M. et al. A high-throughput three-dimensional cell migration assay for toxicity screening with mobile device-based macroscopic image analysis. Sci. Rep. 3, 3000 (2013).
3. Kovács, M., Tóth, J., Hetényi, C., Málnási-Csizmadia, A. & Sellers, J. R. Mechanism of blebbistatin inhibition of myosin II. J. Biol. Chem. 279, 35557-63 (2004).

Scan the QR-codes in the poster for videos of the printed structures and our publication in Scientific Reports!

Acknowledgements/Funding

This work was supported by the National Science Foundation [Small Business Innovation Research Phase I 0945954 and Phase II 1127551]; and the State of Texas Emerging Technology Fund. The authors would like to thank Ahmed Nasser Salameh, Ph.D., Alexes Daquinag, Ph.D., and Mikhail Kolonin, Ph.D., University of Texas Health Science Center at Houston, for donating the 3T3s